



ELK Biotechnology

Histone H2B Rabbit pAb
Catalog Number: ES2526
For research use only.

Overview

Product Name	Histone H2B polyclonal antibody
Host species	Rabbit
Applications	WB,IHC-p,IF,ELISA
Species Cross-Reactivity	Human,Mouse
Recommended dilutions	Western Blot:1/500-1/2000 Immunohistochemistry:1/100-1/300 Immunofluorescence:1/200-1/1000 ELISA:1/5000
Immunogen	Synthesized peptide derived from the N-terminal region of human Histone H2B. at AA range; 10-90
Specificity	Histone H2B Polyclonal Antibody detects endogenous levels of Histone H2B protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	Histone H2B
Gene Name	HIST1H2BH
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	18kDa
Human Gene ID	8345
Human Swiss-Prot Number	Q93079
Alternative Names	HIST1H2BH; H2BFJ; Histone H2B type 1-H; Histone H2B.j; H2B/j