



# ELK Biotechnology

**HSC 70 Rabbit pAb**  
**Catalog Number: ES2570**  
**For research use only.**

## Overview

---

Product Name	HSC 70 polyclonal antibody
Host species	Rabbit
Applications	<b>WB,ELISA</b>
Species Cross-Reactivity	<b>Human,Mouse,Rat</b>
Recommended dilutions	<b>Western Blot:1/500-1/2000</b> <b>ELISA:1/10000</b>
Immunogen	Synthesized peptide derived from the Internal region of human HSC 70.
Specificity	HSC 70 Polyclonal Antibody detects endogenous levels of HSC 70 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	HSC 70
Gene Name	HSPA8
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	<b>72kDa</b>
Human Gene ID	3312
Human Swiss-Prot Number	P11142
Alternative Names	HSPA8; HSC70; HSP73; HSPA10; Heat shock cognate 71 kDa protein; Heat shock 70 kDa protein 8