



ELK Biotechnology

IKK γ Rabbit pAb
Catalog Number: ES2604
For research use only.

Overview

Product Name	IKK γ polyclonal antibody
Host species	Rabbit
Applications	WB,IHC-p,ELISA
Species Cross-Reactivity	Human
Recommended dilutions	Western Blot:1/500-1/2000 Immunohistochemistry:1/100-1/300 ELISA:1/20000
Immunogen	Synthesized peptide derived from human IKK γ around the non-phosphorylation site of S31.
Specificity	IKK γ Polyclonal Antibody detects endogenous levels of IKK γ protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	IKK γ
Gene Name	IKBKG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	48kDa
Human Gene ID	8517
Human Swiss-Prot Number	Q9Y6K9
Alternative Names	IKBKG; FIP3; NEMO; NF-kappa-B essential modulator; NEMO; FIP-3; I κ B kinase-associated protein 1; IKKAP1; Inhibitor of nuclear factor kappa-B kinase subunit gamma; I-kappa-B kinase subunit gamma; IKK-gamma; IKKG; I κ B kinase subunit gamma; NF