



# ELK Biotechnology

**Inhibin  $\beta$ -B Rabbit pAb**  
**Catalog Number: ES2622**  
**For research use only.**

## Overview

Product Name	Inhibin $\beta$ -B polyclonal antibody
Host species	Rabbit
Applications	<b>WB,ELISA</b>
Species Cross-Reactivity	<b>Human,Mouse,Rat</b>
Recommended dilutions	<b>Western Blot:1/500-1/2000</b> <b>ELISA:1/10000</b>
Immunogen	Synthesized peptide derived from the Internal region of human Inhibin $\beta$ -B.
Specificity	Inhibin $\beta$ -B Polyclonal Antibody detects endogenous levels of Inhibin $\beta$ -B protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	Inhibin $\beta$ -B
Gene Name	INHBB
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	<b>45kDa</b>
Human Gene ID	3625
Human Swiss-Prot Number	P09529
Alternative Names	INHBB; Inhibin beta B chain; Activin beta-B chain