

Integrin β5 rabbit pAb

Cat No.:ES2633

For research use only

Overview

Product Name	Integrin β5 rabbit pAb
Host species	Rabbit
Applications	WB;ELISA
Species Cross-Reactivity	Human;Mouse;Rat;Monkey
Recommended dilutions	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not
	yet tested in other applications.
Immunogen	The antiserum was produced against synthesized
	peptide derived from human Integrin beta5. AA
	range:721-770
Specificity	Integrin β5 Polyclonal Antibody detects endogenous
	levels of Integrin β5 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and
	0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Integrin beta-5
Gene Name	ITGB5
Cellular localization	Membrane; Single-pass type I membrane protein.
Purification	The antibody was affinity-purified from rabbit
	antiserum by affinity-chromatography using
	epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	88kD
Human Gene ID	3693
Human Swiss-Prot Number	P18084
Alternative Names	ITGB5; Integrin beta-5
Background	function:Integrin alpha-V/beta-5 is a receptor for
	fibronectin. It recognizes the sequence R-G-D in its
	ligand., similarity: Belongs to the integrin beta chain
	family., similarity: Contains 1 PSI
	domain.,similarity:Contains 1 VWFA
	domain.,subunit:Heterodimer of an alpha and a
	beta subunit. Beta-5 associates with alpha-V.,




+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C







ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C