



# ELK Biotechnology

**Ku-70 Rabbit pAb**  
**Catalog Number: ES2682**  
**For research use only.**

## Overview

<b>Product Name</b>	Ku-70 polyclonal antibody
<b>Host species</b>	Rabbit
<b>Applications</b>	<b>WB,IHC-p,ELISA</b>
<b>Species Cross-Reactivity</b>	<b>Human,Mouse</b>
<b>Recommended dilutions</b>	<b>Western Blot:1/500-1/2000</b> <b>Immunohistochemistry:1/100-1/300</b> <b>ELISA:1/5000</b>
<b>Immunogen</b>	Synthesized peptide derived from human Ku-70 around the non-phosphorylation site of S5.
<b>Specificity</b>	Ku-70 Polyclonal Antibody detects endogenous levels of Ku-70 protein.
<b>Formulation</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
<b>Protein Name</b>	Ku-70
<b>Gene Name</b>	XRCC6
<b>Purification</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Clonality</b>	Polyclonal
<b>Concentration</b>	1 mg/ml
<b>Observed band</b>	<b>70kDa</b>
<b>Human Gene ID</b>	2547
<b>Human Swiss-Prot Number</b>	P12956
<b>Alternative Names</b>	XRCC6; G22P1; X-ray repair cross-complementing protein 6; 5'-deoxyribose-5-phosphate lyase Ku70; 5'-dRP lyase Ku70; 70 kDa subunit of Ku antigen; ATP-dependent DNA helicase 2 subunit 1; ATP-dependent DNA helicase II 70 kDa subunit; CTC box-