

KV1.5 rabbit pAb

Cat No.:ES2688

For research use only

Overview

Product Name	KV1.5 rabbit pAb
Host species	Rabbit
Applications	WB;ELISA;IHC
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:500-2000;IHC-p 1:50-300; ELISA 2000-20000
Immunogen	The antiserum was produced against synthesized
-	peptide derived from human KCNA5. AA
	range:253-302
Specificity	KV1.5 Polyclonal Antibody detects endogenous
	levels of KV1.5 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and
	0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Potassium voltage-gated channel subfamily A
	member 5
Gene Name	KCNA5
Cellular localization	Cell membrane ; Multi-pass membrane protein .
Purification	The antibody was affinity-purified from rabbit
	antiserum by affinity-chromatography using
	epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	68kD
Human Gene ID	3741
Human Swiss-Prot Number	P22460
Alternative Names	KCNA5; Potassium voltage-gated channel subfamily
	A member 5; HPCN1; Voltage-gated potassium
	channel HK2; Voltage-gated potassium channel
	subunit Kv1.5
Background	Potassium channels represent the most complex
	class of voltage-gated ino channels from both
	functional and structural standpoints. Their diverse
	functions include regulating neurotransmitter



+86-27-59760950

ELKbio@ELKbiotech.com


www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C



release, heart rate, insulin secretion, neuronal excitability, epithelial electrolyte transport, smooth muscle contraction, and cell volume. Four sequence-related potassium channel genes - shaker, shaw, shab, and shal - have been identified in Drosophila, and each has been shown to have human homolog(s). This gene encodes a member of the potassium channel, voltage-gated, shaker-related subfamily. This member contains six membrane-spanning domains with a shaker-type repeat in the fourth segment. It belongs to the delayed rectifier class, the function of which could restore the resting membrane potential of beta cells after depolarization and thereby contribute to the regulation of

Western Blot analysis of various cells using KV1.5 Polyclonal Antibody diluted at 1:500



Western blot analysis of KCNA5 Antibody. The lane on the right is blocked with the KCNA5 peptide.



+86-27-59760950


www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C



94KD 66KD 45KD 35KD 26KD Western blot analysis of 293T with KV1.5 Rabbit pAb diluted at 1:2,000.

Immunohistochemical analysis of paraffin-embedded human Gastric adenocarcinoma. 1, Antibody was diluted at 1:200(4° overnight). 2, Tris-EDTA,pH9.0 was used for antigen retrieval. 3,Secondary antibody was diluted at 1:200(room temperature, 45min).





+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C