



ELK Biotechnology

LKB1 Rabbit pAb
Catalog Number: ES2717
For research use only.

Overview

| | |
|---------------------------------|--|
| Product Name | LKB1 polyclonal antibody |
| Host species | Rabbit |
| Applications | WB,ELISA |
| Species Cross-Reactivity | Human,Mouse |
| Recommended dilutions | Western Blot:1/500-1/2000 ELISA:1/5000 |
| Immunogen | Synthesized peptide derived from human LKB1 around the non-phosphorylation site of T189. |
| Specificity | LKB1 Polyclonal Antibody detects endogenous levels of LKB1 protein. |
| Formulation | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Storage | Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles. |
| Protein Name | LKB1 |
| Gene Name | STK11 |
| Purification | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Clonality | Polyclonal |
| Concentration | 1 mg/ml |
| Observed band | 65kDa |
| Human Gene ID | 6794 |
| Human Swiss-Prot Number | Q15831 |
| Alternative Names | STK11; LKB1; PJS; Serine/threonine-protein kinase STK11; Liver kinase B1; LKB1; hLKB1; Renal carcinoma antigen NY-REN-19 |