



ELK Biotechnology

LKB1 Rabbit pAb
Catalog Number: ES2718
For research use only.

Overview

Product Name	LKB1 polyclonal antibody
Host species	Rabbit
Applications	WB,ELISA
Species Cross-Reactivity	Human,Mouse
Recommended dilutions	Western Blot:1/500-1/2000 ELISA:1/20000
Immunogen	Synthesized peptide derived from human LKB1 around the non-phosphorylation site of S334.
Specificity	LKB1 Polyclonal Antibody detects endogenous levels of LKB1 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	LKB1
Gene Name	STK11
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	52kDa
Human Gene ID	6794
Human Swiss-Prot Number	Q15831
Alternative Names	STK11; LKB1; PJS; Serine/threonine-protein kinase STK11; Liver kinase B1; LKB1; hLKB1; Renal carcinoma antigen NY-REN-19