

MEK-1/2 rabbit pAb

Cat No.:ES2772

For research use only

Overview

Product Name MEK-1/2 rabbit pAb

Host species Rabbit

Applications IF;WB;IHC;ELISA

Species Cross-Reactivity Human; Mouse; Rat; Monkey

Recommended dilutions IF: 1:50-200 WB 1:500-2000, ELISA 1:10000-20000

IHC 1:50-300

Immunogen The antiserum was produced against synthesized

peptide derived from human MEK1/2. AA

range:189-238

Specificity MEK-1/2 Polyclonal Antibody detects endogenous

levels of MEK-1/2 protein.

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Storage Store at -20°C. Avoid repeated freeze-thaw cycles.

Protein Name Dual specificity mitogen-activated protein kinase

kinase 1/2

Gene Name MAP2K1/MAP2K2

Cellular localization Cytoplasm, cytoskeleton, microtubule organizing

center, centrosome . Cytoplasm, cytoskeleton, microtubule organizing center, spindle pole body . Cytoplasm . Nucleus . Membrane ; Peripheral

membrane protein. Localizes at centrosomes during

prometaphase, midzone during anaphase and

midbody during telophase/cytokinesis

(PubMed:14737111). Membrane localization is probably regulated by its interaction with KSR1

(PubMed:10409742). .

Purification The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

ClonalityPolyclonalConcentration1 mg/mlObserved band43kD





Human Gene ID Alternative Names


Human Swiss-Prot Number

5604/5605 Q02750/P36507

MAP2K1; MEK1; PRKMK1; Dual specificity mitogen-activated protein kinase kinase 1; MAP kinase kinase 1; MAPKK 1; MKK1; ERK activator kinase 1; MAPK/ERK kinase 1; MEK 1; MAP2K2; MEK2; MKK2; PRKMK2; Dual specificity mitogen-activated protein k

Background


The protein encoded by this gene is a member of the dual specificity protein kinase family, which acts as a mitogen-activated protein (MAP) kinase kinase. MAP kinases, also known as extracellular signal-regulated kinases (ERKs), act as an integration point for multiple biochemical signals. This protein kinase lies upstream of MAP kinases and stimulates the enzymatic activity of MAP kinases upon wide variety of extra- and intracellular signals. As an essential component of MAP kinase signal transduction pathway, this kinase is involved in many cellular processes such as proliferation, differentiation, transcription regulation and development. [provided by RefSeq, Jul 2008],



Immunofluorescence analysis of Hela cell. 1,MEK-1/2 Polyclonal Antibody(red) was diluted at 1:200(4° overnight). CYCS Monoclonal Antibody(4B10)(green) was diluted at 1:200(4° overnight). 2, Goat Anti Rabbit Alexa Fluor 594 Catalog:RS3611 was diluted at 1:








Immunofluorescence analysis of rat-lung tissue.

1,MEK-1/2 Polyclonal Antibody(red) was diluted at

1:200(4°C,overnight). 2, Cy3 labled Secondary antibody
was diluted at 1:300(room temperature, 50min).3, Picture
B: DAPI(blue) 10min. Picture A:Target. Picture B: DAPI.
Picture C: merge of A+B




Immunofluorescence analysis of rat-lung tissue.

1,MEK-1/2 Polyclonal Antibody(red) was diluted at

1:200(4°C,overnight). 2, Cy3 labled Secondary antibody
was diluted at 1:300(room temperature, 50min).3, Picture

B: DAPI(blue) 10min. Picture A:Target. Pictur



Immunofluorescence analysis of mouse-lung tissue.

1,MEK-1/2 Polyclonal Antibody(red) was diluted at

1:200(4°C,overnight). 2, Cy3 labled Secondary antibody
was diluted at 1:300(room temperature, 50min).3, Picture
B: DAPI(blue) 10min. Picture A:Target. Picture B: DAPI.
Picture C: merge of A+B

