



ELK Biotechnology

NOS3 Rabbit pAb
Catalog Number: ES2976
For research use only.

Overview

| | |
|--------------------------|---|
| Product Name | NOS3 polyclonal antibody |
| Host species | Rabbit |
| Applications | WB,IHC-p,IF,ELISA |
| Species Cross-Reactivity | Human,Mouse,Rat |
| Recommended dilutions | Western Blot:1/500-1/2000 Immunohistochemistry:1/100-1/300 Immunofluorescence:1/200-1/1000 ELISA:1/40000 |
| Immunogen | Synthesized peptide derived from human NOS3 around the non-phosphorylation site of S1177. |
| Specificity | NOS3 Polyclonal Antibody detects endogenous levels of NOS3 protein. |
| Formulation | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Storage | Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles. |
| Protein Name | NOS3 |
| Gene Name | NOS3 |
| Purification | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Clonality | Polyclonal |
| Concentration | 1 mg/ml |
| Observed band | 140kDa |
| Human Gene ID | 4846 |
| Human Swiss-Prot Number | P29474 |
| Alternative Names | NOS3; Nitric oxide synthase; endothelial; Constitutive NOS; cNOS; EC-NOS; Endothelial NOS; eNOS; NOS type III; NOSIII |