



ELK Biotechnology

Oct-1 Rabbit pAb
Catalog Number: ES2999
For research use only.

Overview

Product Name	Oct-1 polyclonal antibody
Host species	Rabbit
Applications	WB,IHC-p,ELISA
Species Cross-Reactivity	Human,Mouse,Rat
Recommended dilutions	Western Blot:1/500-1/2000 Immunohistochemistry:1/100-1/300 ELISA:1/20000
Immunogen	Synthesized peptide derived from human 41183 around the non-phosphorylation site of S385.
Specificity	Oct-1 Polyclonal Antibody detects endogenous levels of 40817 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	Oct-1
Gene Name	POU2F1
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	100kDa
Human Gene ID	5451/5452
Human Swiss-Prot Number	P14859/P09086
Alternative Names	POU2F1; OCT1; OTF1; POU domain; class 2, transcription factor 1; NF-A1; Octamer-binding protein 1; Oct-1; Octamer-binding transcription factor 1; OTF-1; POU2F2; OCT2; OTF2; POU domain, class 2, transcription factor 2; Lymphoid-restricted im