

Frizzled-10 Rabbit pAb Catalog Number: ES4010 For research use only.

Overview

Product Name Frizzled-10 polyclonal antibody

Host species Rabbit

Applications WB,IHC-p,ELISA

Species Cross-Reactivity Human, Mouse

Recommended dilutions Western Blot:1/500-1/2000

Immunohistochemistry:1/100-1/300

ELISA:1/20000

Immunogen Synthesized peptide derived from the N-terminal region of human Frizzled-

10.

Specificity Frizzled-10 Polyclonal Antibody detects endogenous levels of Frizzled-10

protein.

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid

repeated freeze / thaw cycles.

Protein Name Frizzled-10

Gene Name FZD10

Purification The antibody was affinity-purified from rabbit antiserum by affinity-

chromatography using epitope-specific immunogen.

Clonality Polyclonal

Concentration 1 mg/ml

Observed band 65kDa

Human Gene ID 11211

Human Swiss-Prot Number Q9ULW2

Alternative Names FZD10; Frizzled-10; Fz-10; hFz10; FzE7; CD350