



ELK Biotechnology

GDF-9 Rabbit pAb
Catalog Number: ES4031
For research use only.

Overview

Product Name	GDF-9 polyclonal antibody
Host species	Rabbit
Applications	WB,ELISA
Species Cross-Reactivity	Human,Mouse,Rat
Recommended dilutions	Western Blot:1/500-1/2000 ELISA:1/20000
Immunogen	Synthesized peptide derived from the C-terminal region of human GDF-9.
Specificity	GDF-9 Polyclonal Antibody detects endogenous levels of GDF-9 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	GDF-9
Gene Name	GDF9
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	51kDa
Human Gene ID	2661
Human Swiss-Prot Number	O60383
Alternative Names	GDF9; Growth/differentiation factor 9; GDF-9