



ELK Biotechnology

BMP-8A Rabbit pAb
Catalog Number: ES4036
For research use only.

Overview

| | |
|---------------------------------|---|
| Product Name | BMP-8A polyclonal antibody |
| Host species | Rabbit |
| Applications | WB,IHC-p,ELISA |
| Species Cross-Reactivity | Human,Mouse,Rat |
| Recommended dilutions | Western Blot:1/500-1/2000 Immunohistochemistry:1/100-1/300 ELISA:1/20000 |
| Immunogen | Synthesized peptide derived from the Internal region of human BMP-8A. |
| Specificity | BMP-8A Polyclonal Antibody detects endogenous levels of BMP-8A protein. |
| Formulation | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Storage | Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles. |
| Protein Name | BMP-8A |
| Gene Name | BMP8A |
| Purification | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Clonality | Polyclonal |
| Concentration | 1 mg/ml |
| Observed band | 44kDa |
| Human Gene ID | 353500 |
| Human Swiss-Prot Number | Q7Z5Y6 |
| Alternative Names | BMP8A; Bone morphogenetic protein 8A; BMP-8A |