



ELK Biotechnology

Lyn Rabbit pAb
Catalog Number: ES4043
For research use only.

Overview

Product Name	Lyn polyclonal antibody
Host species	Rabbit
Applications	WB,IHC-p,ELISA
Species Cross-Reactivity	Human,Mouse,Rat
Recommended dilutions	Western Blot:1/500-1/2000 Immunohistochemistry:1/100-1/300 ELISA:1/20000
Immunogen	Synthesized peptide derived from human Lyn around the non-phosphorylation site of Y192.
Specificity	Lyn Polyclonal Antibody detects endogenous levels of Lyn protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	Lyn
Gene Name	LYN
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	58kDa
Human Gene ID	4067
Human Swiss-Prot Number	P07948
Alternative Names	LYN; JTK8; Tyrosine-protein kinase Lyn; Lck/Yes-related novel protein tyrosine kinase; V-yes-1 Yamaguchi sarcoma viral related oncogene homolog; p53Lyn; p56Lyn