



ELK Biotechnology

NFκB-p65 Rabbit pAb
Catalog Number: ES4069
For research use only.

Overview

| | |
|---------------------------------|--|
| Product Name | NFκB-p65 polyclonal antibody |
| Host species | Rabbit |
| Applications | WB,IHC-p,ELISA |
| Species Cross-Reactivity | Human,Mouse,Rat |
| Recommended dilutions | Western Blot:1/500-1/2000 Immunohistochemistry:1/100-1/300 ELISA:1/20000 |
| Immunogen | Synthesized peptide derived from human NFκB-p65 around the non-acetylation site of K218. |
| Specificity | NFκB-p65 Polyclonal Antibody detects endogenous levels of NFκB-p65 protein. |
| Formulation | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Storage | Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles. |
| Protein Name | NFκB-p65 |
| Gene Name | RELA |
| Purification | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Clonality | Polyclonal |
| Concentration | 1 mg/ml |
| Observed band | 60kDa |
| Human Gene ID | 5970 |
| Human Swiss-Prot Number | Q04206 |
| Alternative Names | RELA; NFKB3; Transcription factor p65; Nuclear factor NF-kappa-B p65 subunit; Nuclear factor of kappa light polypeptide gene enhancer in B-cells 3 |