



## ELK Biotechnology

### MMP-9 Rabbit pAb

Catalog Number: ES4087

For research use only.

#### Overview

---

Product Name	MMP-9 polyclonal antibody
Host species	Rabbit
Applications	WB,IHC-p,ELISA
Species Cross-Reactivity	Human
Recommended dilutions	Western Blot:1/500-1/2000 Immunohistochemistry:1/100-1/300 ELISA:1/20000
Immunogen	Synthesized peptide derived from the Internal region of human MMP-9.
Specificity	MMP-9 Polyclonal Antibody detects endogenous levels of MMP-9 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	MMP-9
Gene Name	MMP9
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	80kDa
Human Gene ID	4318
Human Swiss-Prot Number	P14780
Alternative Names	MMP9; CLG4B; Matrix metalloproteinase-9; MMP-9; 92 kDa gelatinase; 92 kDa type IV collagenase; Gelatinase B; GELB