



# ELK Biotechnology

**Histone H3 Rabbit pAb**  
**Catalog Number: ES4163**  
**For research use only.**

## Overview

---

<b>Product Name</b>	Histone H3 polyclonal antibody
<b>Host species</b>	Rabbit
<b>Applications</b>	<b>WB,ELISA</b>
<b>Species Cross-Reactivity</b>	<b>Human,Mouse,Rat</b>
<b>Recommended dilutions</b>	<b>Western Blot:1/500-1/2000</b> <b>ELISA:1/20000</b>
<b>Immunogen</b>	Synthesized peptide derived from human Histone H3 around the non-methylation site of K10.
<b>Specificity</b>	Histone H3 Polyclonal Antibody detects endogenous levels of Histone H3 protein only when non-methylated at K10.
<b>Formulation</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
<b>Protein Name</b>	Histone H3
<b>Gene Name</b>	HIST1H3A/HIST1H3/HIST1H3C/HIST1H3D/HIST1H3E/HIST1H3F/HIST1H3G/HIST1H3H/HIST1H3I/HIST1H3J/HIST2H3A/HIST2H3C/HIST2H3D/H3F3A/H3F3B/H3F3C
<b>Purification</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Clonality</b>	Polyclonal
<b>Concentration</b>	1 mg/ml
<b>Observed band</b>	<b>kDa</b>
<b>Human Gene ID</b>	8350/8351/8352/8353/8354/8355/8356/8357/8358/8968/126961/333932/653604/3020/3021/440093
<b>Human Swiss-Prot Number</b>	P68431/Q71DI3/P84243/Q6NXT2
<b>Alternative Names</b>	HIST1H3A; H3FA; HIST1H3B; H3FL; HIST1H3C; H3FC; HIST1H3D; H3FB; HIST1H3E; H3FD; HIST1H3F; H3FI; HIST1H3G; H3FH; HIST1H3H; H3FK; HIST1H3I; H3FF; HIST1H3J; H3FJ; Histone H3.1; Histone H3.1; Histone H3.2; Histone H3/m; Histone H3/o; H3F3A; H3.3A; H3F3; PP781;