



ELK Biotechnology

Integrin β 4 Rabbit pAb
Catalog Number: ES4198
For research use only.

Overview

Product Name	Integrin β 4 polyclonal antibody
Host species	Rabbit
Applications	WB,ELISA
Species Cross-Reactivity	Human,Mouse,Rat
Recommended dilutions	Western Blot:1/500-1/2000 ELISA:1/20000
Immunogen	Synthesized peptide derived from the Internal region of human Integrin β 4.
Specificity	Integrin β 4 Polyclonal Antibody detects endogenous levels of Integrin β 4 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	Integrin β 4
Gene Name	ITGB4
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	200kDa
Human Gene ID	3691
Human Swiss-Prot Number	P16144
Alternative Names	ITGB4; Integrin beta-4; GP150; CD104