



ELK Biotechnology

HDAC1 Rabbit pAb
Catalog Number: ES4225
For research use only.

Overview

Product Name	HDAC1 polyclonal antibody
Host species	Rabbit
Applications	WB,IHC-p,ELISA
Species Cross-Reactivity	Human,Mouse,Rat
Recommended dilutions	Western Blot:1/500-1/2000 Immunohistochemistry:1/100-1/300 ELISA:1/20000
Immunogen	Synthesized peptide derived from the Internal region of human HDAC1.
Specificity	HDAC1 Polyclonal Antibody detects endogenous levels of HDAC1 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	HDAC1
Gene Name	HDAC1
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	55kDa
Human Gene ID	3065
Human Swiss-Prot Number	Q13547
Alternative Names	HDAC1; RPD3L1; Histone deacetylase 1; HD1