



# ELK Biotechnology

**SRY Rabbit pAb**  
**Catalog Number: ES4227**  
**For research use only.**

## Overview

---

Product Name	SRY polyclonal antibody
Host species	Rabbit
Applications	<b>WB,IHC-p,ELISA</b>
Species Cross-Reactivity	<b>Human</b>
Recommended dilutions	<b>Western Blot:1/500-1/2000</b> <b>Immunohistochemistry:1/100-1/300</b> <b>ELISA:1/20000</b>
Immunogen	Synthesized peptide derived from the Internal region of human SRY.
Specificity	SRY Polyclonal Antibody detects endogenous levels of SRY protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	SRY
Gene Name	SRY
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	<b>25kDa</b>
Human Gene ID	6736
Human Swiss-Prot Number	Q05066
Alternative Names	SRY; TDF; Sex-determining region Y protein; Testis-determining factor