



ELK Biotechnology

SREBP-1 Rabbit pAb
Catalog Number: ES4238
For research use only.

Overview

Product Name	SREBP-1 polyclonal antibody
Host species	Rabbit
Applications	WB,ELISA
Species Cross-Reactivity	Human,Mouse,Rat
Recommended dilutions	Western Blot:1/500-1/2000 ELISA:1/20000
Immunogen	Synthesized peptide derived from the Internal region of human SREBP-1.
Specificity	SREBP-1 Polyclonal Antibody detects endogenous levels of SREBP-1 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	SREBP-1
Gene Name	SREBF1
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	120kDa
Human Gene ID	6720
Human Swiss-Prot Number	P36956
Alternative Names	SREBF1; BHLHD1; SREBP1; Sterol regulatory element-binding protein 1; SREBP-1; Class D basic helix-loop-helix protein 1; bHLHd1; Sterol regulatory element-binding transcription factor 1