



ELK Biotechnology

IGFBP3 Rabbit pAb
Catalog Number: ES4248
For research use only.

Overview

Product Name	IGFBP3 polyclonal antibody
Host species	Rabbit
Applications	WB,IHC-p,ELISA
Species Cross-Reactivity	Human,Mouse,Rat
Recommended dilutions	Western Blot:1/500-1/2000 Immunohistochemistry:1/100-1/300 ELISA:1/20000
Immunogen	Synthesized peptide derived from the Internal region of human IGFBP3.
Specificity	IGFBP3 Polyclonal Antibody detects endogenous levels of IGFBP3 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	IGFBP3
Gene Name	IGFBP3
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	30kDa
Human Gene ID	3486
Human Swiss-Prot Number	P17936
Alternative Names	IGFBP3; IBP3; Insulin-like growth factor-binding protein 3; IBP-3; IGF-binding protein 3; IGFBP-3