



## ELK Biotechnology

**Cortactin Rabbit pAb**  
**Catalog Number: ES4269**  
**For research use only.**

### Overview

---

|                          |   |
|--------------------------|---|
| Product Name             | Cortactin polyclonal antibody   |
| Host species             | Rabbit  |
| Applications             | <b>WB,ELISA</b>   |
| Species Cross-Reactivity | <b>Human,Mouse,Rat</b>  |
| Recommended dilutions    | <b>Western Blot:1/500-1/2000</b><br><b>ELISA:1/10000</b>  |
| Immunogen                | Synthesized peptide derived from the Internal region of human Cortactin.  |
| Specificity              | Cortactin Polyclonal Antibody detects endogenous levels of Cortactin protein.   |
| Formulation              | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.   |
| Storage                  | Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.               |
| Protein Name             | Cortactin   |
| Gene Name                | CTTN  |
| Purification             | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Clonality                | Polyclonal  |
| Concentration            | 1 mg/ml   |
| Observed band            | <b>62kDa</b>  |
| Human Gene ID            | 2017  |
| Human Swiss-Prot Number  | Q14247  |
| Alternative Names        | CTTN; EMS1; Src substrate cortactin; Amplaxin; Oncogene EMS1  |