



ELK Biotechnology

PSA Rabbit pAb
Catalog Number: ES4281
For research use only.

Overview

Product Name	PSA polyclonal antibody
Host species	Rabbit
Applications	WB,IHC-p,ELISA
Species Cross-Reactivity	Human
Recommended dilutions	Western Blot:1/500-1/2000 Immunohistochemistry:1/100-1/300 ELISA:1/10000
Immunogen	Synthesized peptide derived from the Internal region of human PSA.
Specificity	PSA Polyclonal Antibody detects endogenous levels of PSA protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	PSA
Gene Name	KLK3
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	29kDa
Human Gene ID	354
Human Swiss-Prot Number	P07288
Alternative Names	KLK3; APS; Prostate-specific antigen; PSA; Gamma-seminoprotein; Semin; Kallikrein-3; P-30 antigen; Semenogelase