



# ELK Biotechnology

**RGS1 Rabbit pAb**  
**Catalog Number: ES4289**  
**For research use only.**

## Overview

---

Product Name	RGS1 polyclonal antibody
Host species	Rabbit
Applications	<b>WB,ELISA</b>
Species Cross-Reactivity	<b>Human,Mouse,Rat</b>
Recommended dilutions	<b>Western Blot:1/500-1/2000</b> <b>ELISA:1/10000</b>
Immunogen	Synthesized peptide derived from the Internal region of human RGS1.
Specificity	RGS1 Polyclonal Antibody detects endogenous levels of RGS1 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	RGS1
Gene Name	RGS1
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	<b>23kDa</b>
Human Gene ID	5996
Human Swiss-Prot Number	Q08116
Alternative Names	RGS1; 1R20; BL34; IER1; Regulator of G-protein signaling 1; RGS1; B-cell activation protein BL34; Early response protein 1R20