



NOXIN rabbit pAb

Cat No.:ES4294

For research use only

Overview


| | |
|---------------------------------|--|
| Product Name | NOXIN rabbit pAb |
| Host species | Rabbit |
| Applications | WB;ELISA |
| Species Cross-Reactivity | Human;Rat;Mouse; |
| Recommended dilutions | Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications. |
| Immunogen | The antiserum was produced against synthesized peptide derived from the Internal region of human NOXIN. AA range:111-160 |
| Specificity | NOXIN Polyclonal Antibody detects endogenous levels of NOXIN protein. |
| Formulation | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Storage | Store at -20°C. Avoid repeated freeze-thaw cycles. |
| Protein Name | Nitric oxide-inducible gene protein |
| Gene Name | NOXIN |
| Cellular localization | Cytoplasm . Nucleus . Accumulates in the nucleus in response to stress. . |
| Purification | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Clonality | Polyclonal |
| Concentration | 1 mg/ml |
| Observed band | 85kD |
| Human Gene ID | 220042 |
| Human Swiss-Prot Number | Q8IXT1 |
| Alternative Names | NOXIN; C11orf82; Nitric oxide-inducible gene protein |
| Background | function:Acts as an anti-apoptotic factor and its absence increases cell death under normal and stress conditions. Can induce cell cycle arrest in the G1 or early S phase and this activity is independent of TP53/p53.,subcellular location:Accumulates in |





ELK Biotechnology

the nucleus in response to stress.,



Western Blot analysis of 823, 293T, AD293, Hela, HepG2 cells using NOXIN Polyclonal Antibody. Antibody was diluted at 1:2000. Secondary antibody(catalog#:RS0002) was diluted at 1:20000



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road,Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C