



ELK Biotechnology

NGF Rabbit pAb
Catalog Number: ES4295
For research use only.

Overview

Product Name	NGF polyclonal antibody
Host species	Rabbit
Applications	WB,IHC-p,ELISA
Species Cross-Reactivity	Human
Recommended dilutions	Western Blot:1/500-1/2000 ELISA 1/10000-1/20000 Immunohistochemistry:1/50-1/300
Immunogen	Synthesized peptide derived from the Internal region of human NGF.
Specificity	NGF Polyclonal Antibody detects endogenous levels of NGF protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	NGF
Gene Name	NGF
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	27kDa
Human Gene ID	4803
Human Swiss-Prot Number	P01138
Alternative Names	NGF; NGFB; Beta-nerve growth factor; Beta-NGF