



ELK Biotechnology

ICAM-2 Rabbit pAb
Catalog Number: ES4331
For research use only.

Overview

Product Name	ICAM-2 polyclonal antibody
Host species	Rabbit
Applications	WB,ELISA
Species Cross-Reactivity	Human
Recommended dilutions	Western Blot:1/500-1/2000 ELISA:1/10000
Immunogen	Synthesized peptide derived from the Internal region of human ICAM-2.
Specificity	ICAM-2 Polyclonal Antibody detects endogenous levels of ICAM-2 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	ICAM-2
Gene Name	ICAM2
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	30kDa
Human Gene ID	3384
Human Swiss-Prot Number	P13598
Alternative Names	ICAM2; Intercellular adhesion molecule 2; ICAM-2; CD102