



ELK Biotechnology

NKp30 Rabbit pAb
Catalog Number: ES4361
For research use only.

Overview

Product Name	NKp30 polyclonal antibody
Host species	Rabbit
Applications	WB,ELISA
Species Cross-Reactivity	Human,Rat
Recommended dilutions	Western Blot:1/500-1/2000 ELISA:1/10000
Immunogen	Synthesized peptide derived from the Internal region of human NKp30.
Specificity	NKp30 Polyclonal Antibody detects endogenous levels of NKp30 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	NKp30
Gene Name	NCR3
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	21kDa
Human Gene ID	259197
Human Swiss-Prot Number	O14931
Alternative Names	NCR3; 1C7; LY117; Natural cytotoxicity triggering receptor 3; Activating natural killer receptor p30; Natural killer cell p30-related protein; NK-p30; NKp30; CD337