

TANK rabbit pAb

Cat No.: ES4403

For research use only

Overview

Product Name TANK rabbit pAb

Host species Rabbit

Applications WB;IHC;IF;ELISA Species Cross-Reactivity Human;Mouse;Rat

Recommended dilutions Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.

Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications. The antiserum was produced against synthesize.

Immunogen The antiserum was produced against synthesized

peptide derived from human TANK. AA

range:171-220

Specificity TANK Polyclonal Antibody detects endogenous levels

of TANK protein.

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Storage Store at -20°C. Avoid repeated freeze-thaw cycles.

Protein Name TRAF family member-associated NF-kappa-B

activator

Gene Name TANK

Cellular localization Cytoplasm.

Purification The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal
Concentration 1 mg/ml
Observed band 48kD
Human Gene ID 10010
Human Swiss-Prot Number Q92844

Alternative Names TANK; ITRAF; TRAF2; TRAF family

member-associated NF-kappa-B activator;

TRAF-interacting protein; I-TRAF

Background The TRAF (tumor necrosis factor receptor-associated

factor) family of proteins associate with and




+86-27-59760950 ELKbio@ELKbiotech.com www.




transduce signals from members of the tumor necrosis factor receptor superfamily. The protein encoded by this gene is found in the cytoplasm and can bind to TRAF1, TRAF2, or TRAF3, thereby inhibiting TRAF function by sequestering the TRAFs in a latent state in the cytoplasm. For example, the protein encoded by this gene can block TRAF2 binding to LMP1, the Epstein-Barr virus transforming protein, and inhibit LMP1-mediated NF-kappa-B activation. Three alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Nov 2010],

Western Blot analysis of HT-29 cells using TANK Polyclonal Antibody




Western blot analysis of lysates from HT29 cells, using I-TRAF Antibody. The lane on the right is blocked with the synthesized peptide.



+86-27-59760950







Immunohistochemical analysis of paraffin-embedded human tonsil. 1, Antibody was diluted at 1:200(4° overnight). 2, Tris-EDTA,pH9.0 was used for antigen retrieval. 3,Secondary antibody was diluted at 1:200(room temperature, 30min).



+86-27-59760950