



ELK Biotechnology

p57 Rabbit pAb
Catalog Number: ES4454
For research use only.

Overview

Product Name	p57 polyclonal antibody
Host species	Rabbit
Applications	WB,IF,ELISA
Species Cross-Reactivity	Human,Mouse
Recommended dilutions	Western Blot:1/500-1/2000 Immunofluorescence:1/200-1/1000 ELISA:1/5000
Immunogen	Synthesized peptide derived from human p57 around the non-phosphorylation site of T310.
Specificity	p57 Polyclonal Antibody detects endogenous levels of p57 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	p57
Gene Name	CDKN1C
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	32kDa
Human Gene ID	1028
Human Swiss-Prot Number	P49918
Alternative Names	CDKN1C; KIP2; Cyclin-dependent kinase inhibitor 1C; Cyclin-dependent kinase inhibitor p57; p57Kip2