

## Epac rabbit pAb

## Cat No.:ES4469

For research use only

## Overview

Product Name	Epac rabbit pAb
Host species	Rabbit
Applications	WB;ELISA
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not
	yet tested in other applications.
Immunogen	The antiserum was produced against synthesized
	peptide derived from human RAPGEF3. AA
	range:712-761
Specificity	Epac Polyclonal Antibody detects endogenous levels
	of Epac protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and
	0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Rap guanine nucleotide exchange factor 3
Gene Name	RAPGEF3
Cellular localization	Endomembrane system .
Purification	The antibody was affinity-purified from rabbit
	antiserum by affinity-chromatography using
	epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	100kD
Human Gene ID	10411
Human Swiss-Prot Number	O95398
Alternative Names	RAPGEF3; CGEF1; EPAC; EPAC1; Rap guanine
	nucleotide exchange factor 3; Exchange factor
	directly activated by cAMP 1; Exchange protein
	directly activated by cAMP 1; EPAC 1; Rap1
	guanine-nucleotide-exchange factor directly
	activated by cAMP;
Background	caution:It is uncertain whether Met-1 or Met-43 is
	the initiator.,domain:The DEP domain is involved in



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C



membrane localization independent from regulation by cAMP.,function:Guanine nucleotide exchange factor (GEF) for RAP1A and RAP2A small GTPases that is activated by binding cAMP.,similarity:Contains 1 cyclic nucleotide-binding domain.,similarity:Contains 1 DEP domain.,similarity:Contains 1 N-terminal Ras-GEF domain.,similarity:Contains 1 Ras-GEF domain.,tissue specificity:Widely expressed with highest levels in adult kidney, heart, thyroid and brain, and fetal kidney.,

Western blot analysis of the lysates from HeLa cells using RAPGEF3 antibody.





+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C