



ELK Biotechnology

Cdx2 Rabbit pAb
Catalog Number: ES4474
For research use only.

Overview

| | |
|---------------------------------|---|
| Product Name | Cdx2 polyclonal antibody |
| Host species | Rabbit |
| Applications | WB,ELISA |
| Species Cross-Reactivity | Human,Mouse,Rat |
| Recommended dilutions | Western Blot:1/500-1/2000 ELISA:1/5000 |
| Immunogen | Synthesized peptide derived from the Internal region of human Cdx2. at AA range; 20-100 |
| Specificity | Cdx2 Polyclonal Antibody detects endogenous levels of Cdx2 protein. |
| Formulation | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Storage | Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles. |
| Protein Name | Cdx2 |
| Gene Name | CDX2 |
| Purification | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Clonality | Polyclonal |
| Concentration | 1 mg/ml |
| Observed band | 34kDa |
| Human Gene ID | 1045 |
| Human Swiss-Prot Number | Q99626 |
| Alternative Names | CDX2; CDX3; Homeobox protein CDX-2; CDX-3; Caudal-type homeobox protein 2 |