



ELK Biotechnology

ZnT-9 Rabbit pAb
Catalog Number: ES4480
For research use only.

Overview

Product Name	ZnT-9 polyclonal antibody
Host species	Rabbit
Applications	WB,ELISA
Species Cross-Reactivity	Human,Mouse,Rat
Recommended dilutions	Western Blot:1/500-1/2000 ELISA:1/40000
Immunogen	Synthesized peptide derived from the Internal region of human ZnT-9. at AA range; 140-220
Specificity	ZnT-9 Polyclonal Antibody detects endogenous levels of ZnT-9 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	ZnT-9
Gene Name	SLC30A9
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	60kDa
Human Gene ID	10463
Human Swiss-Prot Number	Q9Y6R2
Alternative Names	