



ELK Biotechnology

C/EBP β Rabbit pAb
Catalog Number: ES4491
For research use only.

Overview

Product Name	C/EBP β polyclonal antibody
Host species	Rabbit
Applications	WB,IHC-p,ELISA
Species Cross-Reactivity	Human,Mouse
Recommended dilutions	Western Blot:1/500-1/2000 Immunohistochemistry:1/100-1/300 ELISA:1/5000
Immunogen	Synthesized peptide derived from the C-terminal region of human C/EBP β . at AA range; 200-280
Specificity	C/EBP β Polyclonal Antibody detects endogenous levels of C/EBP β protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	C/EBP β
Gene Name	CEBPB
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	49kDa
Human Gene ID	1051
Human Swiss-Prot Number	P17676
Alternative Names	CEBPB; LAP; TCF5; PP9092; CCAAT/enhancer-binding protein beta; C/EBP beta; Liver activator protein; Nuclear factor NF-IL6; Transcription factor 5; TCF-5