



ELK Biotechnology

C/EBP ϵ Rabbit pAb
Catalog Number: ES4501
For research use only.

Overview

Product Name	C/EBP ϵ polyclonal antibody
Host species	Rabbit
Applications	WB,IHC-p,ELISA
Species Cross-Reactivity	Human,Mouse,Rat
Recommended dilutions	Western Blot:1/500-1/2000 Immunohistochemistry:1/100-1/300 ELISA:1/10000
Immunogen	Synthesized peptide derived from the C-terminal region of human C/EBP ϵ . at AA range; 160-240
Specificity	C/EBP ϵ Polyclonal Antibody detects endogenous levels of C/EBP ϵ protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	C/EBP ϵ
Gene Name	CEBPE
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	26kDa
Human Gene ID	1053
Human Swiss-Prot Number	Q15744
Alternative Names	CEBPE; CCAAT/enhancer-binding protein epsilon; C/EBP epsilon