



ELK Biotechnology

NTE Rabbit pAb
Catalog Number: ES4570
For research use only.

Overview

Product Name	NTE polyclonal antibody
Host species	Rabbit
Applications	WB,ELISA
Species Cross-Reactivity	Human,Mouse
Recommended dilutions	Western Blot:1/500-1/2000 ELISA:1/20000
Immunogen	Synthesized peptide derived from the N-terminal region of human NTE. at AA range; 80-160
Specificity	NTE Polyclonal Antibody detects endogenous levels of NTE protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	NTE
Gene Name	PNPLA6
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	150kDa
Human Gene ID	10908
Human Swiss-Prot Number	Q8IY17
Alternative Names	PNPLA6; NTE; Neuropathy target esterase; Patatin-like phospholipase domain-containing protein 6