



ELK Biotechnology

MIA2 Rabbit pAb

Catalog Number: ES4664

For research use only.

Overview

Product Name	MIA2 polyclonal antibody
Host species	Rabbit
Applications	WB,ELISA
Species Cross-Reactivity	Human
Recommended dilutions	Western Blot:1/500-1/2000 ELISA:1/20000
Immunogen	Synthesized peptide derived from the Internal region of human MIA2. at AA range; 100-180
Specificity	MIA2 Polyclonal Antibody detects endogenous levels of MIA2 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	MIA2
Gene Name	MIA2
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	62kDa
Human Gene ID	117153
Human Swiss-Prot Number	Q96PC5
Alternative Names	MIA2; Melanoma inhibitory activity protein 2