

PC-PLD4 rabbit pAb

Cat No.:ES4701

For research use only

Overview

Product Name	PC-PLD4 rabbit pAb
Host species	Rabbit
Applications	WB;IHC
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:500-2000;IHC-p 1:50-300
Immunogen	The antiserum was produced against synthesized
-	peptide derived from human PLD4. AA
	range:457-506
Specificity	PC-PLD4 Polyclonal Antibody detects endogenous
. ,	levels of PC-PLD4 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and
	0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Phospholipase D4
Gene Name	PLD4
Cellular localization	Endoplasmic reticulum membrane ; Single-pass type
	II membrane protein . Golgi apparatus, trans-Golgi
	network membrane ; Single-pass type II membrane
	protein . Nucleus . Early endosome . Cytoplasmic
	vesicle, phagosome . Activation of microglia induces
	translocation of PLD4 from the nucleus to the
	phagosomes
Purification	The antibody was affinity-purified from rabbit
	antiserum by affinity-chromatography using
	epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	50kD
Human Gene ID	122618
Human Swiss-Prot Number	Q96BZ4
Alternative Names	PLD4; C14orf175; Phospholipase D4; PLD 4; Choline
	phosphatase 4; Phosphatidylcholine-hydrolyzing
	phospholipase D4



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C



Background

Jurkat

PC-PLD4

-- 117

-- 85

-- 48 -- 34

-- 26

-- 19 (kD)

PLD4

138-100-70--55-

40--35--25--

15-

PLD4 --

(kD)

117-

85-

48-

34-

26-

19-

catalytic activity:A phosphatidylcholine + H(2)O = choline + a phosphatidate.,similarity:Belongs to the phospholipase D family.,similarity:Contains 2 PLD phosphodiesterase domains.,

Western Blot analysis of Jurkat cells using PC-PLD4 Polyclonal Antibody

Western blot analysis of lysates from Jurkat cells, using PLD4 Antibody. The lane on the right is blocked with the synthesized peptide.


Western blot analysis of the lysates from K562 cells using PLD4 antibody.



www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C





Immunohistochemical analysis of paraffin-embedded human Squamous cell carcinoma of lung. 1, Antibody was diluted at 1:200(4° overnight). 2, Tris-EDTA,pH9.0 was used for antigen retrieval. 3,Secondary antibody was diluted at 1:200(room temperature, 45min).



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C