

CNG-1 rabbit pAb

Cat No.:ES4717

For research use only

Overview

| Product Name | CNG-1 rabbit pAb |
|------------------------------|---|
| Host species | Rabbit |
| Applications | WB;ELISA |
| Species Cross-Reactivity | Human; Mouse; Rat |
| Recommended dilutions | Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not |
| | yet tested in other applications. |
| Immunogen | The antiserum was produced against synthesized |
| | peptide derived from human CNGA1. AA |
| | range:401-450 |
| Specificity | CNG-1 Polyclonal Antibody detects endogenous |
| | levels of CNG-1 protein. |
| Formulation | Liquid in PBS containing 50% glycerol, 0.5% BSA and |
| | 0.02% sodium azide. |
| Storage | Store at -20°C. Avoid repeated freeze-thaw cycles. |
| Protein Name | cGMP-gated cation channel alpha-1 |
| Gene Name | CNGA1 |
| Cellular localization | Cell membrane ; Multi-pass membrane protein . |
| Purification | The antibody was affinity-purified from rabbit |
| | antiserum by affinity-chromatography using |
| | epitope-specific immunogen. |
| Clonality | Polyclonal |
| Concentration | 1 mg/ml |
| Observed band | 80kD |
| Human Gene ID | 1259 |
| Human Swiss-Prot Number | P29973 |
| Alternative Names | CNGA1; CNCG; CNCG1; cGMP-gated cation channel |
| | alpha-1; Cyclic nucleotide-gated cation channel 1; |
| | Cyclic nucleotide-gated channel alpha-1; CNG |
| | channel alpha-1; CNG-1; CNG1; Cyclic |
| | nucleotide-gated channel; photoreceptor; Rod |
| | photoreceptor cG |
| Background | The protein encoded by this gene is involved in |
| | phototransduction. Along with another protein, the |



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C



HepG2

CNG-1

138---

70---55---40---35---

15----

CNGA1

(kD)

117-

85-

48-

34-

26-

19-

HepG2 K562

--100

--75 --50 --37

--25 --20

--15 (kd)

CNGA1

+86-27-59760950

encoded protein forms a cGMP-gated cation channel in the plasma membrane, allowing depolarization of rod photoreceptors. This represents the last step in the phototransduction pathway. Defects in this gene are a cause of retinitis pigmentosa autosomal recessive (ARRP) disease. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Dec 2008],

Western Blot analysis of HepG2 cells using CNG-1 Polyclonal Antibody

Western blot analysis of CNGA1 Antibody. The lane on the right is blocked with the CNGA1 peptide.

Western blot analysis of the lysates from HepG2 cells using CNGA1 antibody.

www.elkbiotech.com



23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C

ELKbio@ELKbiotech.com





+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C