



Gab 2/4 rabbit pAb

Cat No.:ES4743

For research use only

Overview


Product Name	Gab 2/4 rabbit pAb
Host species	Rabbit
Applications	WB;IHC;IF;ELISA
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Immunogen	The antiserum was produced against synthesized peptide derived from human GAB2. AA range:21-70
Specificity	Gab 2/4 Polyclonal Antibody detects endogenous levels of Gab 2/4 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	GRB2-associated-binding protein 2/4
Gene Name	GAB2/GAB4
Cellular localization	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	62kD
Human Gene ID	128954
Human Swiss-Prot Number	Q2WGN9/Q9UQC2
Alternative Names	GAB4; GRB2-associated-binding protein 4; GRB2-associated binder 2-like; GAB2-like; GRB2-associated binder 4; GRB2-associated-binding protein 2-like; Growth factor receptor bound protein 2-associated protein 4; GAB2; KIAA0571; GRB2-associate
Background	similarity:Belongs to the GAB family.,similarity:Contains 1 PH domain.,





ELK Biotechnology

Immunohistochemistry analysis of paraffin-embedded human tesophagus, using GAB4 Antibody. The picture on the right is blocked with the synthesized peptide.



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road,Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C