



ELK Biotechnology

GPR175 Rabbit pAb
Catalog Number: ES4762
For research use only.

Overview

Product Name	GPR175 polyclonal antibody
Host species	Rabbit
Applications	IF,ELISA
Species Cross-Reactivity	Human
Recommended dilutions	Immunofluorescence:1/200-1/1000 ELISA:1/5000
Immunogen	Synthesized peptide derived from the C-terminal region of human GPR175. at AA range; 240-320
Specificity	GPR175 Polyclonal Antibody detects endogenous levels of GPR175 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	GPR175
Gene Name	TPRA1
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	kDa
Human Gene ID	131601
Human Swiss-Prot Number	Q86W33
Alternative Names	TPRA1; GPR175; PP6566; Transmembrane protein adipocyte-associated 1; Integral membrane protein GPR175