



ELK Biotechnology

Histone H1 Rabbit pAb
Catalog Number: ES4766
For research use only.

Overview

Product Name	Histone H1 polyclonal antibody
Host species	Rabbit
Applications	IHC-p,IF,ELISA
Species Cross-Reactivity	Human
Recommended dilutions	Immunohistochemistry:1/100-1/300 Immunofluorescence:1/200-1/1000 ELISA:1/10000
Immunogen	Synthesized peptide derived from human Histone H1 around the non-acetylation site of K25. at AA range; 100-180
Specificity	Histone H1 Polyclonal Antibody detects endogenous levels of Histone H1 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	Histone H1
Gene Name	H1FOO
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	kDa
Human Gene ID	132243
Human Swiss-Prot Number	Q8IZA3
Alternative Names	H1FOO; H1OO; OSH1; Histone H1oo; Oocyte-specific histone H1; Oocyte-specific linker histone H1; osH1