



ELK Biotechnology

GPR151 Rabbit pAb
Catalog Number: ES4774
For research use only.

Overview

Product Name	GPR151 polyclonal antibody
Host species	Rabbit
Applications	WB,IF,ELISA
Species Cross-Reactivity	Human
Recommended dilutions	Western Blot:1/500-1/2000 Immunofluorescence:1/200-1/1000 ELISA:1/20000
Immunogen	Synthesized peptide derived from the C-terminal region of human GPR151. at AA range; 350-430
Specificity	GPR151 Polyclonal Antibody detects endogenous levels of GPR151 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	GPR151
Gene Name	GPR151
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	47kDa
Human Gene ID	134391
Human Swiss-Prot Number	Q8TDV0
Alternative Names	GPR151; PGR7; Probable G-protein coupled receptor 151; G-protein coupled receptor PGR7; GPCR-2037