



# ELK Biotechnology

## CREB-1 Rabbit pAb Catalog Number: ES4800 For research use only.

### Overview

---

Product Name	CREB-1 polyclonal antibody
Host species	Rabbit
Applications	WB,IHC-p,ELISA
Species Cross-Reactivity	Human,Mouse,Rat
Recommended dilutions	<b>Western Blot:1/500-1/2000</b> <b>Immunohistochemistry:1/100-1/300</b> <b>ELISA:1/20000</b>
Immunogen	Synthesized peptide derived from human CREB-1 around the non-phosphorylation site of T100.
Specificity	CREB-1 Polyclonal Antibody detects endogenous levels of CREB-1 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	CREB-1
Gene Name	CREB1
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	<b>36kDa</b>
Human Gene ID	1385
Human Swiss-Prot Number	P16220
Alternative Names	CREB1; Cyclic AMP-responsive element-binding protein 1; CREB-1; cAMP-responsive element-binding protein 1