

## CD116 Rabbit pAb Catalog Number: ES4849 For research use only.

## Overview

Product Name CD116 polyclonal antibody

Host species Rabbit

Applications WB,ELISA

Species Cross-Reactivity Human

Recommended dilutions Western Blot:1/500-1/2000

ELISA:1/20000

**Immunogen** Synthesized peptide derived from the Internal region of human CD116.

**Specificity** CD116 Polyclonal Antibody detects endogenous levels of CD116 protein.

**Formulation** Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid

repeated freeze / thaw cycles.

Protein Name CD116

Gene Name CSF2RA

Purification The antibody was affinity-purified from rabbit antiserum by affinity-

chromatography using epitope-specific immunogen.

**Clonality** Polyclonal

Concentration 1 mg/ml

Observed band 46kDa

Human Gene ID 1438

Human Swiss-Prot Number P15509

Alternative Names CSF2R4; CSF2R; CSF2RY; Granulocyte-macrophage colony-stimulating

factor receptor subunit alpha; GM-CSF-R-alpha; GMCSFR-alpha; GMR-

alpha; CDw116; CD antigen CD116