

CMTM4 rabbit pAb

Cat No.: ES4867

For research use only

Overview

Product Name CMTM4 rabbit pAb

Host species Rabbit

Applications WB;IHC;IF;ELISA **Species Cross-Reactivity** Human;Mouse

Recommended dilutions Western Blot: 1/500 - 1/2000. IHC-p: 1:100-300

ELISA: 1/20000. IF 1:100-300 Not yet tested in other

applications.

Immunogen The antiserum was produced against synthesized

peptide derived from human CKLF4. AA

range:161-210

Specificity CMTM4 Polyclonal Antibody detects endogenous

levels of CMTM4 protein.

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Store at -20°C. Avoid repeated freeze-thaw cycles.

Protein Name CKLF-like MARVEL transmembrane

domain-containing protein 4

Gene Name CMTM4

Cellular localizationMembrane ; Multi-pass membrane protein .PurificationThe antibody was affinity-purified from rabbitanticorum by affinity shromatography using

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal
Concentration 1 mg/ml
Observed band 26kD
Human Gene ID 146223
Human Swiss-Prot Number Q8IZR5

Alternative Names CMTM4; CKLFSF4; CKLF-like MARVEL

transmembrane domain-containing protein 4; Chemokine-like factor superfamily member 4

Background This gene belongs to the chemokine-like factor gene


superfamily, a novel family that is similar to the chemokine and the transmembrane 4 superfamilies




+86-27-59760950 ELKbio@ELKbiotech.com www




of signaling molecules. This gene is one of several chemokine-like factor genes located in a cluster on chromosome 16. Alternatively spliced transcript variants encoding different isoforms have been identified. [provided by RefSeq, Jul 2008],



Immunofluorescence analysis of A549. 1,primary
Antibody(red) was diluted at 1:200(4°C overnight). 2, Goat
Anti Rabbit IgG (H&L) - Alexa Fluor 594 Secondary
antibody was diluted at 1:1000(room temperature,
50min).3, Picture B: DAPI(blue) 10min.



Western Blot analysis of HELA cells using CMTM4 Polyclonal Antibody



Western blot analysis of CKLF4 Antibody. The lane on the right is blocked with the CKLF4 peptide.




+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com





Immunohistochemistryt analysis of paraffin-embedded human breast carcinoma, using CKLF4 Antibody. The lane on the right is blocked with the CKLF4 peptide.



+86-27-59760950