



ELK Biotechnology

GDF-7 Rabbit pAb
Catalog Number: ES4907
For research use only.

Overview

Product Name	GDF-7 polyclonal antibody
Host species	Rabbit
Applications	WB,ELISA
Species Cross-Reactivity	Human
Recommended dilutions	Western Blot:1/500-1/2000 ELISA:1/40000
Immunogen	Synthesized peptide derived from the Internal region of human GDF-7. at AA range; 160-240
Specificity	GDF-7 Polyclonal Antibody detects endogenous levels of GDF-7 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	GDF-7
Gene Name	GDF7
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	47kDa
Human Gene ID	151449
Human Swiss-Prot Number	Q7Z4P5
Alternative Names	GDF7; Growth/differentiation factor 7; GDF-7